

SKINTYTE

SKINTYTE

For
finer, more
youthful
looking skin.



Before

After 1 treatment



Photos courtesy of Dr. Gloria Lopez



Untreated right hand

Left hand after 1 treatment

Photos courtesy of Dr. Gloria Lopez



SKINTYTE

SkinTyte™ offers a non-ablative approach to firmer skin.

What is SkinTyte?

The SkinTyte treatment is based on the principle of selective thermolysis and utilizes Sciton's BBL™ pulsed light technology to deeply heat your skin while cooling and protecting the entire treated surface. SkinTyte can be provided quickly and comfortably with no downtime. The end result is a more youthful appearance.

How does SkinTyte work?

The SkinTyte procedure utilizes a broadband of light energy to deeply heat dermal collagen. This process promotes partial denaturation and contraction of collagen. The heat also initiates the body's natural healing process which creates a renewed collagen foundation leading to increased skin firmness. Powerful cooling maintains the outer surface of your skin at a cool temperature before, during, and after each pulse providing the patient a comfortable and safe procedure.



Photos courtesy of Dr. Jason Pozner

Before

After 3 treatments

What is the treatment like?

The SkinTyte light energy is delivered in a sequence of rapid, gentle pulses. There is no need for a topical anesthetic, however, your physician may choose to use one on more sensitive areas. A series of pulses will be repeated over the treatment area to ensure best results. While the face and neck can typically be treated within 30 minutes, actual time will depend on the area treated.

What should I expect after treatment?

You may experience some mild redness in the treated area immediately after and this should return to normal within minutes. Since the SkinTyte procedure is non-invasive there is no downtime so you can resume regular activities immediately following treatment. You can also apply makeup and sunscreen after each session.

What results should I expect?

The SkinTyte procedure is not a substitute for a face-lift and results may vary from patient to patient. The SkinTyte procedure is recommended for patients who desire moderate, noticeable skin improvement without having to undergo surgery. Full results often take up to four months as the collagen foundation of your skin is rebuilt and strengthened.

What areas can be treated?

SkinTyte treatments are safe for all skin types and can be performed on any area of the body where an improvement in your skin firmness is desired. Popular areas include the face, neck, abdomen, and arms.

You should consult with your physician regarding the treatment options that are best for you.



Before



Immediately after

Photos courtesy of Dr. Gloria Lopez

How many treatments will I need?

A series of 2-3 treatments spaced a month apart is typical of most patients. Your physician will personalize a treatment plan based on the level of improvement you desire. Your physician may also combine SkinTyte with other popular aesthetic procedures to further personalize your treatment.

What aftercare do I need?

No special skincare regimen is required following treatment. However, you should always protect your skin from direct sun exposure and use a good sunscreen. Your physician can also recommend an at home skin care regimen if desired.

©2007 Sciton, Inc.

Printed in Australia.

BBL and SkinTyte are trademarks of Sciton, Inc. Sciton is a registered mark. 2600-026-00 Rev. A